



DIFFERENTIAL CELL COUNTERS

Each key has the name and picture of a cell type above it, representing particular cell groups. The key label strips are easily removed to allow customisation, if required.

Simple and intuitive differential cell totalising counters to simplify blood cell counts.

**A bell rings to indicate when total of 100 cells has been counted.
Percentage of each cell type is then read directly in window above the relevant key.
All counts are simply zeroed by counter-clockwise rotation of the end knobs.**

8 PLACE and 5 PLACE versions.

Strong, reliable keys
Solid and sturdy construction
UK Manufacturer

Order Code: [WDC 035](#)

8 PLACE: ([as illustration](#))

to count Basophils, Eosinophils, Myelophils, Juven., Stab., Seg., Lymphocytes, Monocytes

Order Code: [WDC 025](#)

5 PLACE :

to count Basophils, Eosinophils, Seg., Lymphocytes, Monocytes

